
kliniderm® foam sacral

INSTRUCTIONS FOR USE

KLINION

Distributor: Medico B.V. | Alexander Flemingstraat 2 | NL-2261 MA Oud-Beerland | The Netherlands

Advanced Medical Sciences Ltd. | 1remier Park 33 Road Ore, Winford Industrial Estate

NOINITK

ENGLISH

Product Description
Kliniderm Foam Sacral dressing is a sterile wound dressing, consisting of a pink waterproof polyurethane bacterial barrier film, coated with perforated pressure sensitive acrylic adhesive with a centrally located hydrophilic absorbent polyurethane foam pad. The wound contact layer of the dressing is a perforated film (coated on both sides with a pressure sensitive acrylic adhesive) which protects the wound bed by preventing the formation of granulation tissue into the pores of the foam, thus reducing trauma on dressing removal. The perforations in the wound contact layer allow the uptake of exudate and prevent excess fluid causing maceration to healthy skin. Kliniderm Foam Sacral dressing maintains a moist environment to assist in the healing process of moderate to heavily exuding wounds.

Indications

Kliniderm Foam Sacral dressing is indicated for use on moderate to heavily exuding chronic and acute wounds. The dressing may be used throughout the healing process on the following wounds; pressure ulcers, venous leg ulcers, diabetic ulcers, lacerations, abrasions, post operative surgical wounds, superficial and partial thickness burns and cavity wounds, as a secondary dressing. Kliniderm Foam Sacral dressings are suitable for use under compression bandaging.

Directions for Use

1. Irrigate the wound site in accordance with standard protocols if necessary and dry the surrounding skin before applying the dressing.
2. Select a size of Kliniderm Foam Sacral dressing making sure the foam pad size is slightly larger than the wound.
3. The foam pad size should extend well beyond the open wound [at least 2cm] in all directions.
4. Remove the release liners from the wound dressing.
5. Centre the dressing on the wound, pink side up and apply it gently to the wound site. Press the dressing on the surrounding skin by gentle application of pressure.

Dressing Change and Removal

1. Kliniderm Foam Sacral dressing can remain in situ up to 7 days, dependent on patient condition and the level of exudate. Change the dressing when significant signs of exudate are visible on the dressing. Initially, it may be necessary to observe the wound or dressing frequently and change the dressing whenever good wound care practice dictates.
2. Remove the dressing by peeling off from a corner using gentle pressure.
3. Gently remove the dressing from the wound bed and dispose of according to local procedures and guidelines.
4. If required, cleanse the wound site with a suitable wound cleanser prior to application of a new dressing.

Precautions

Do not use on individuals with a known sensitivity to polyurethane films or foams or acrylic adhesives. Do not reuse in whole or in part, as it may compromise sterility and/or the performance of the dressing. Do not use with oxidising solutions such as hypochlorite or hydrogen peroxide, as these can break down the absorbent polyurethane component of the dressing. Kliniderm Foam Sacral dressings must not be ingested and must be kept away from children and animals. For use under the guidance of a healthcare professional. For external use only.

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