
kliniderm® silicone wound contact layer

INSTRUCTIONS FOR USE

KLINION®

Medeco B.V. | Alexander Flemingstraat 2 | NL-3261 MA Oud-Beijerland | The Netherlands

KLINION®

ENGLISH

Product Description

Silicone wound contact layer is a transparent, flexible wound dressing. It is non-adherent to wound site, minimizes trauma to the wound and pain to the patient during dressing changes. The silicone wound contact layer consists of 3 layers: a PET film, a polyurethane net that is one-side coated with soft silicon and a release film. The middle layer is the functional component and the coated surface should contact the wound. The porous structure of the wound contact layer allows exudate to pass into an outer absorbent dressing. The wound contact layer prevents the outer dressing from sticking to the wound and ensures atraumatic dressing changes. The wound contact layer seals around the wound edges, preventing the exudate from leaking onto the surrounding skin, thus minimizing the risk of maceration. High transparency allows wound inspection during application and wear.

Indications

Silicone contact layer can be used with a secondary absorbent dressing, indicated for a wide range of exuding wounds, such as diabetic foot ulcer, pressure ulcer, leg venous and arterial ulcer, first and second degree burns, surgical wounds and traumatic wounds. It can also be used as wound protective layer on non-exuding wounds and on areas with fragile skin. The product is intended for a professional use only. Product may remain in place for up to 14 days which allows cost effective and undisturbed wound healing.

Contra-indications

- 1) Do not use on Third-degree burns.
- 2) Do not use on the patients with a known hypersensitivity to the product itself or to its components.
- 3) Not for surgical implantation.

Directions for use

First step: Clean the wound area in accordance with normal procedures. Ensure the peri-wound skin is dry. Second step: Products shall be kept clean and sterile before use, do not use if package is damaged. Select a size of silicone contact layer enough to overlap the wound edges by at least 1cm. large wounds may require more overlap. Third step: Remove the release film and gently apply the wound contact layer directly on the wound site; peel off the dressing carrier and smooth down the dressings edges to ensure good adhesion. Fourth step: Apply outer absorbent dressing and a retention dressing if needed.

Frequency of change

The silicone wound contact layer should be changed according to clinical condition of the wound. Silicone wound contact layer may be left to the wound up to 14 days. If the condition of the wound deteriorates, discontinue use.

Precautions

1. The wound should be inspected for signs of infection according to clinical practice. Infected wounds should be regularly inspected and only be treated by a healthcare professional.
2. If the silicone contact layer is used on Epidermolysis Bullosa patients, employ extra surveillance at dressing changes.
3. If the silicone contact layer is used on burns treated with meshed grafts: avoid placing unnecessary pressure upon the dressing.
4. If the silicone contact layer is used after facial resurfacing: avoid placing pressure upon the dressing. Lift and reposition the dressing at least every second day.
5. If the silicone contact layer is used on bleeding wounds or wounds with high viscosity exudate, the dressing should be covered with a moist absorbent dressing pad.
6. When the silicone contact layer is used for the fixation of skin grafts, the dressing should not be changed before the fifth day post application.
7. Do not reuse. Reusing may cause cross-contamination and serious infection
8. If the silicone contact layer is used on burns treated with meshed grafts or after facial resurfacing, imprints may occur if the product is not used properly.

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