# L-Mesitran® OINTMENT Hydro-active antibacterial wound ointment

NETERLANDS

L-Mesitran Zalf bevat: 48% medicinale honing, hypoallergene lanoline (Medilan®), zonnebloemolie, levertraanolie, Calendula Officinalis, Aloe Barbadensis, vitamine C & E en zinkoxide.

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ENGLISH

L-Mesitran Ointment contains: 48% medical grade honey, medical grade hypoallergenic lanolin (Medilan®), sunflower oil, cod liver oil, Calendula Officinalis, Aloe Barbadensis, vitamin C & E and zinc oxide.

When L-Mesitran is in contact with the wound, fluid is drawn from the surrounding tissues, this creates a moist wound healing environment. This stimulates the wound healing process by facilitating autolysis of necrotic

and devitalised material, aids wound bed preparation and promotes epithelization.

L-Mesitran:

- has antibacterial properties.
  - dilutes gradually and can be removed easily as per local wound cleansing protocol.
  - prevents the secondary dressing adhering to the wound bed.
- Due to the use of natural ingredients, each batch may vary in colour (yellow - brown), viscosity and may contain visible sugar crystals. These variances will not affect the efficacy and use of the product.

## Indications

Superficial and acute wounds e.g. cuts, abrasions and donor sites. Superficial and partial thickness burns (1st & 2nd degree). Chronic wounds e.g. pressure ulcers, venous, arterial and diabetic ulcers. Fungating wounds (to help deodorise and debride). Colonised acute wounds and (postoperative) surgical wounds.

## Contraindications and constraints

No contraindications are known to date. Do not use on patients who are sensitive to the product or any of its components.

## Preventive measures and comments

1. Always consult a doctor for wound management advice.
2. Do not use new products that are damaged or expired.
3. Inspect and clean the wound according to local protocol.
4. Due to autolytic debridement the wound may appear to increase in size after the first dressing changes. This is normal and occurs as the devitalised tissue is removed, preparing the wound bed for healing.
5. In the event of intolerance to L-Mesitran or one of its components, remove it using appropriate techniques and clean the area carefully.
6. L-Mesitran is for single patient use, the product may only be re-used with the same patient.
7. Always re-seal the product after use.
8. After opening use the product within three months.
9. Dispose of the product as per local protocol.

## Instructions for use

A. Tube: 1. Remove the lid, clean the mouth of the tube and the top of the lid with an alcohol swab. 2. Push the top of the lid into the foil seal that covers the nozzle, do not squeeze the tube at this time. B. Jar: remove the seal and lid. **Application to the wound:** 1. Use a sterile spatula to distribute a thin layer of L-Mesitran either directly to the wound area and the surrounding tissue or on to sterile gauze first. 2. Cover and secure with the appropriate secondary dressing depending on the amount of wound fluids. 3. Inspect the wound and the surrounding tissue regularly, as per local protocols.

## Frequency of use

This depends upon the amount of exudate and the needs of the wound. Re-apply the product when it has dissolved (± once every 24-48 hours).

## Storage precautions

Store at room temperature between 5°C and 25°C (41°-77°F). Avoid direct sunlight and moisture. L-Mesitran is a natural product and it may become liquid under higher temperatures. If this is expected to be a problem place the product in a refrigerator prior to use.

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