



Soft Fine® Plus Pen-Nadeln



KLINION Soft Fine® Plus
0,23 mm (32G) x 4 mm
100 + 10 Stück
PZN 10179649

KLINION Soft Fine® Plus
0,23 mm (32G) x 6 mm
100 + 10 Stück
PZN 11340242

KLINION Soft Fine® Plus
0,23 mm (32G) x 8 mm
100 + 10 Stück
PZN 11340259

KLINION Soft Fine® Plus
0,23 mm (32G) x 5 mm
100 + 10 Stück
PZN 09076312

KLINION Soft Fine® Plus
0,25 mm (31G) x 6 mm
100 + 10 Stück
PZN 09166245

KLINION Soft Fine® Plus
0,25 mm (31G) x 8 mm
100 + 10 Stück
PZN 09166251

KLINION Soft Fine® Plus
0,33 mm (29G) x 12 mm
100 + 10 Stück
PZN 09166274

KLINION® Soft Fine® Plus Pen-Nadeln sind für alle Insulinpens mit Feingewinde geeignet:

Novo Nordisk
NovoPen
Echo



Novo Nordisk
NovoPen 5



sanofi-aventis
TactiPen



sanofi-aventis
ClikStar



Berlin-Chemie
BerliPen
junior



Berlin-Chemie
BerliPen
Precision



Berlin-Chemie
BerliPen
areo 3



Lilly
HumaPen
Savvio



Lilly
HumaPen
Luxura HD



Owen-Mumford
AutoPen 24



Haselmeier
DiaPen



sanofi-aventis - JuniorStar

Novo Nordisk - NovoPen Junior, NovoPen 4

Berlin-Chemie - BerliPen areo, BerliPen areo 2, BerlinPen 301 und 302

Lilly - HumaPen Memoir, HumaPen Luxura

Owen-Mumford - AutoPen classic

Ypsomed - Ypsopen

Pendiq und Smart Pendiq 2.0

eu-medical GmbH • Bärensteiner Str. 27-29 • 01277 Dresden
info@eu-medical.de • www.eu-medical.de • Gebührenfreie Rufnummer: 0800 - 5 88 80 56

Stand: September 2020