

# Kliniderm Fiber



## Hydroactive fibre wound dressing

### Product description

Kliniderm Fiber is a soft, highly absorbent, malleable, gel-forming wound dressing. Thanks to its high (vertical) absorption capacity, wound exudate is absorbed in the dressing and converted to gel. This conversion provides a moist wound environment, which stimulates healing and autolytic debridement. The dressing can be easily removed without damage to the tissue.



Kliniderm Fiber

### Indications

Kliniderm Fiber can be used on moderately to strongly exuding wounds, both chronic and acute (see classification diagram on rear). It is also suitable for stemming small bleeds in surface wounds. Kliniderm Fiber wound dressing can be used together with compression bandages.


### Use of Kliniderm Fiber

This strongly absorbent Kliniderm Fiber wound dressing creates a gel as soon as it comes into contact with wound exudate, and takes on the shape of the wound bed, flexibly and comfortably. Kliniderm Fiber can be left on the wound up to 7 days, depending on the wound condition.


### Absorption capacity



### Tensile strength wet (in use)




### Transverse absorption



# Hydroactive fibre wound dressing

## Classification diagram



### Reduces transverse absorption

Kliniderm Fiber has concentrated, vertical fluid absorption capacity, thus preventing maceration of the wound edges. Spread to the surrounding skin is therefore minimal; certainly when compared with standard alginate wound dressings, for example. Kliniderm Fiber is also excellent for use in compression therapy.



Vertical absorption

### Super strong

Kliniderm Fiber is even three times stronger than hydrofibres in use (during changes). The wound dressing is therefore simple to remove in one piece. Debris is removed without damaging the tissue, upon removing the dressing. The aim is to create a clean wound bed, after all. No residue is left in the wound.

### Great absorption capacity

Kliniderm Fiber wound dressing absorbs a great deal of exudate and creates an ideal, moist wound environment that stimulates wound healing. The dressing is particularly suitable for moderately to strongly exuding wounds.



Great absorption capacity

## ORDERING INFORMATION

REF no.	Dimension	Contents of packaging:	Contents of cover box
174520	5 x 5 cm	10 pcs	10 x 10 x 1 pc
174521	10 x 10 cm	10 pcs	10 x 10 x 1 pc
174522	15 x 15 cm	5 pcs	10 x 5 x 1 pc
174530	42 x 2 cm	5 pcs	10 x 5 x 1 pc