

klinipress gauze & non-woven compresses

When treating wounds, especially chronic wounds, moist wound treatment using special wound dressings is generally the first choice. In some cases however, the treater will opt for dry wound treatment using traditional wound dressings, such as gauze and non-woven compresses. When opting for a traditional product to dress the wound the dressing must be changed daily, sometimes even more often. This ensures the compress does not stick to the wound.

Traditional compresses can also be used for cleaning wounds. For quick cleaning, the compress is moistened with a wound rinsers (NaCl, for example) at room temperature and the wound is dabbed clean. In the case of stubborn wound tissue, the moist compress is applied to the wound for 15 minutes in order to soften it first.

Klinion® gauze and non-woven compresses are available in various sizes both non-sterile and sterile. There are 4 lines available:

High quality line with high absorbency available both in sterile and non-sterile: 69 REF

- Gauze 17 threads and 40s 100% cotton yarn with folded edges
- NW 70% viscose and 30% polyester, 30g/m²

Tenderline both in sterile and non-sterile: 36 REF

- Gauze 13 threads and 40s 100% cotton yarn without folded edges
- NW 50% viscose and 50% polyester, 30g/m²

Hygiene packing line: 8 REF

- Gauze 17 threads and 40s 100% cotton yarn without folded edges
- NW 70% viscose and 30% polyester, 30g/m²

OR gauze line: 17 REF

- Blue Xray thread
- Gauze 17 threads and 36s/40s 100% cotton yarn with folded edges
- Double blister packing to guarantee sterility to highest degree
- EO sterilized
- Soft and low lint
- With counting card



features & benefits

- Low lint
- High absorbency
- Very good range available both sterile and non-sterile
- Double blister packing guarantee sterility to highest degree on OR gauze line
- Also hygiene packing available for more easy and hygienic way of using the products
- Cost effective solution

product details

Intended use	This product is for treating wounds, especially chronic wounds, moist wound treatment using special wound dressings is generally the first choice
Components / materials	100% cotton (gauze). Non-woven dressings 70% viscose/30% polyester or 50% viscose/50% polyester
PVC / phthalates	No
Latex	No
CE classification	1, 1S, 2a
Sterile product	Not sterilized or EO sterilized
Special labeling	NA
Storage unopened	5 years
Storing opened	5 years, not applicable for sterile products
Disposal of the product	Follow local regulations for hazardous waste
Disposal of inner packaging	Follow local regulations for packing material waste
Disposal of outer packaging	Follow local regulations for packing material waste

high quality gauze

ref nr.	dimensions	layers	pcs in blister	pcs in box	pcs in carton
gauze non-sterile					
4111025	5x5cm	8		100	6000
4111030	7,5x7,5cm	8		100	2000
4111098	10x10cm	8		100	2000
40512130	10x20cm	8		100	3000
4111028	5x5cm	12		100	1500
4111029	7,5x7,5cm	12		100	2000
4111101	10x10cm	12		100	1000
4111124	10x20cm	12		50	500
40512135	10x10cm	16		100	3000
40512140	10x20cm	24		50	1000
40512145	20x20cm	32		50	300
gauze sterile					
4111009	5x5cm	8	1	100	6000
4111010	5x5cm	8	2	100	3000
4111011	5x5cm	8	3	150	4500
4111015	5x5cm	8	5	150	3000
4111040	7,5x7,5cm	8	1	100	2000
4111041	7,5x7,5cm	8	2	100	2000
4111042	7,5x7,5cm	8	3	150	3000
4111043	7,5x7,5cm	8	5	150	2250
4111002	8,5x8,5cm	8	1	100	1600
4111082	10x10cm	8	1	100	1600
4111080	10x10cm	8	2	100	2000
4111079	10x10cm	8	3	150	3000
4111078	10x10cm	8	5	150	1500
4111221w	20x20cm	8	2	60	960
4111235w	20x20cm	8	5	50	800
4111017	5x5cm	12	1	100	2000
4111018	5x5cm	12	2	100	2000
4111019	5x5cm	12	3	90	1080
4111020	5x5cm	12	5	100	1200
4111021w	5x5cm	12	10	100	1200
4111022	7,5x7,5cm	12	5	100	4000
4111085	10x10cm	12	1	100	1600
4111086	10x10cm	12	2	100	1600
4111087	10x10cm	12	3	105	1260
4111088	10x10cm	12	5	100	1200
4111089-1	10x10cm	12	10	100	800
4111112	10x20cm	12	1	50	600
4111113	10x20cm	12	2	80	1280
4111114	10x20cm	12	3	90	1440
40512160	20x25cm	12	2	60	240
40512150	10x10cm	16	5	125	1000
4111223w	20x20cm	16	1	30	480
4111255w	20x25cm	16	1	25	400
4111256-1	20x30cm	16	1	20	160
40512155	20x20cm	32	5	30	180

high quality non-woven

ref nr.	dimensions	layers	pcs in blister	pcs in box	pcs in carton
non-woven non-sterile					
4175004	5x5cm	4		100	5000
4175024	7,5x7,5cm	4		100	5000
4175034	10x10cm	4		100	5000
4175043	10x20cm	4		100	2000
4175104	5x5cm	8		100	5000
4175134	10x10cm	8		100	5000
non-woven sterile					
4175001	5x5cm	4	1	100	2000
4175002	5x5cm	4	2	100	2000
4175003	5x5cm	4	5	100	3000
4175008	7,5x7,5cm	4	2	100	2000
4175009	7,5x7,5cm	4	5	100	2000
4175030	10x10cm	4	1	100	1600
4175031	10x10cm	4	2	100	1600
4175033-1	10x10cm	4	5	100	1600
4175040	10x20cm	4	1	50	600
4175041	10x20cm	8	2	50	600
4175101	5x5cm	8	1	100	200
4175102	5x5cm	8	2	100	200
4175110	5x10cm	8	1	16	384
4175130	10x10cm	8	1	100	1600
4175131	10x10cm	8	2	100	1600
4175140	10x20cm	8	1	50	600
4175141	10x20cm	8	2	50	600

tenderline gauze

ref nr.	dimensions	layers	pcs in blister	pcs in box	pcs in carton
gauze non-sterile					
40512080	5x5cm	8		100	24000
40512081	7,5x7,5cm	8		100	12000
40512082	10x10cm	8		100	9000
40512083	10x20cm	8		100	4500
gauze sterile					
40512090	5x5cm	8	2	100	4800
40512100	5x5cm	8	5	100	5250
40512091	7,5x7,5cm	8	2	150	3000
40512101	7,5x7,5cm	8	5	150	3000
40512092	10x10cm	8	2	100	1600
40512102	10x10cm	8	5	100	2000
40512093	10x20cm	8	2	150	1200
40512103	10x20cm	8	5	150	1200
40512110	5x5cm	12	2	100	4800
40512120	5x5cm	12	5	100	5250
40512111	7,5x7,5cm	12	2	100	3000
40512121	7,5x7,5cm	12	5	150	3000
40512112	10x10cm	12	2	150	1600
40512122	10x10cm	12	5	60	2000
40512113	10x20cm	12	2	50	1200
40512123	10x20cm	12	5	100	1200

hygiene packing

ref nr.	dimensions	layers	pcs in blister	pcs in box	pcs in carton
gauze non-sterile					
40512002	5x5cm	8		120	10800
40512003	7,5x7,5cm	8		120	7200
40512004	10x10cm	8		120	7200
40512005	10x20cm	12		50	1000
non-woven non-sterile					
40512022	5x5cm	4		150	13500
40512023	7,5x7,5cm	4		150	9000
40512024	10x10cm	4		150	9000
40512025	10x20cm	4		75	1500

tenderline non-woven

ref nr.	dimensions	layers	pcs in blister	pcs in box	pcs in carton
non-woven non-sterile					
40512050	5x5cm	4		100	24000
40512051	7,5x7,5cm	4		100	12000
40512052	10x10cm	4		100	9000
40512053	10x20cm	4		100	4500
40512030	5x5cm	4		200	5000
40512031	7,5x7,5cm	4		200	5000
40512032	10x10cm	4		200	5000
40512033	10x20cm	4		200	5000
non-woven sterile					
40512060	5x5cm	4	2	160	4800
40512070	5x5cm	4	5	175	5250
40512061	7,5x7,5cm	4	2	200	3000
40512071	7,5x7,5cm	4	5	200	3000
40512062	10x10cm	4	2	200	1600
40512070	10x10cm	4	5	250	2000
40512063	10x20cm	4	2	200	1200
40512073	10x20cm	4	5	200	1200

or gauzes

ref nr.	dimensions	layers	pcs in blister	pcs in box	pcs in carton
gauze sterile					
115001	5x5cm	12	10	100	800
115007	10x10cm	12	10	100	900
115006	10x10cm	12	20	200	1200
115005	10x10cm	12	5	100	900
115021	10x20cm	12	10	100	600
115020	10x20cm	12	5	100	600
115023	10x20cm	24	5	50	300
115032	20x20cm	32	5	25	150
115047	20x25cm	16	2	30	180
115046	20x25cm	16	5	40	240
115048	20x25cm	32	2	16	96
115062	45x45cm	4	5	40	320
115064	45x45cm*	4	5	40	320
115065	45x70cm*	4	2	20	80
115060	45x70cm	4	1	20	80
115066	45x70cm	4	3	21	84
115063	45x70cm	4	5	25	100

*with loop